MISSOURI DEPARTMENT OF HEALTH STATE PUBLIC HEALTH LABORATORY DATAMASTER MAINTENANCE REPORT

BREATH ALCOHOL PROGRAM

HW-9 2001

Complete this report in duplicate at the time of the regular monthly preventive maintenance check, and when is repaired. Send copy to Department of Health; retain original in department file. DATAMASTER SN 204089 Inventory# 127214 LOCATION OF INSTRUMENT ISTREET AND CITY) Stoddard County Sheriff's Department, 101 Court Street, Bloomfield, MO 63825 CHECKLIST: Place a check () to the left of each item if found to be satisfactory or if operating within established in observed values where determined.) Unchecked items must be corrected before using instrument. DIAGNOSTIC CHECK (PRINTOUT ATTACHED)			
Stoddard County Sheriff's Department, 101 Court Street, Bloomfield, MO 63825 CHECKLIST: Place a check (<) to the left of each item if found to be satisfactory or if operating within established in observed values where determined.) Unchecked items must be corrected before using instrument.			
in observed values where determined.) Unchecked items must be corrected before using instrument.			
of the control was the corrected before using instrument.	od limits. (W		
1 1: ((((((((((((((
☑ COMPUTER ☑ DETECTOR			
[7]	☑ FILTERS		
☑ HEATERS SAMPLE CHAMBER 50 °C ☑ QUARTZ STANDARD			
✓ FLOW DETECTOR ✓ CALIBRATION			
☑ PUMP HIGH SPEED ☑ PRINTER			
☑ INDICATOR LIGHTS			
☑ TIME AND DATE			
☑ SIMULATOR TEMPERATURE (34 °C ± 0.2°C) 34.0			
☑ CALIBRATION CHECK -			
Run three tests using a standard solution. All three tests must be within \pm 5% of the standard value and spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTAIRED RECIRCULATION PLIMP)	must have CHED) (USI		
RECIRCULATION PUMP) Ø 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE	must have . CHED) (USI		
RECIRCULATION PUMP) Ø 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)	CHED) (US		
RECIRCULATION PUMP) Ø 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .098	CHED) (US		
RECIRCULATION PUMP) © 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 0.098 TEST 2 0.098 TEST 3 0.099	OHED) (USE		
RECIRCULATION PUMP) O. 100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .098 PERFORM R.F.I. TEST (PRINTOUT ATTACHED) NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS)	OHED) (USE		
RECIRCULATION PUMP) O .100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE O .040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .09 PERFORM R.F.I. TEST (PRINTOUT ATTACHED) NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH 1 FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS) EFUSALS 2 (004) 4 (.0509) 7 (.1014) 1 (.1519) 3 (Over .19) st any new parts and describe any alteration or modification that was made to restore the instrument to	OHED) (USE		
RECIRCULATION PUMP) ② 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE ① 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) EST 1 .098 TEST 2 .098 TEST 3 .09 PERFORM R.F.I. TEST (PRINTOUT ATTACHED) ② NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS) EFUSALS 2 (004) 4 (.0509) 7 (.1014) 1 (.1519) 3 (Over .19) st any new parts and describe any alteration or modification that was made to restore the instrument to	OHED) (USE		
RECIRCULATION PUMP) O .100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE O .040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .09 PERFORM R.F.I. TEST (PRINTOUT ATTACHED) NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS) EFUSALS 2 (004) 4 (.0509) 7 (.1014) 1 (.1519) 3 (Over .19) ist any new parts and describe any alteration or modification that was made to restore the instrument to	OHED) (USE		
RECIRCULATION PUMP) ② 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE ① 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE ② 0.040% STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1	OHED) (USE		
RECIRCULATION PUMP) ② 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE □ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1	OHED) (USE		
RECIRCULATION PUMP) ② 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE □ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .09 ✓ PERFORM R.F.I. TEST (PRINTOUT ATTACHED) ✓ NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS) EFUSALS 2 (004) 4 (.0509) 7 (.1014) 1 (.1519) 3 (Over .19) Ist any new parts and describe any alteration or modification that was made to restore the instrument to operate and within established limits (use other side if necessary) Guth Laboratories Inc., Lot Number 08340, Expiration Date 10/15/2009	OHED) (USE		
RECIRCULATION PUMP) ② 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE □ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) TEST 1 .098 TEST 2 .098 TEST 3 .098 ☑ PERFORM R.F.I. TEST (PRINTOUT ATTACHED) ☑ NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS) EFUSALS 2 (004) 4 (.0509) 7 (.1014) 1 (.1519) 3 (Over.19) ist any new parts and describe any alteration or modification that was made to restore the instrument to operate and within established limits (use other side if necessary) Guth Laboratories Inc., Lot Number 08340, Expiration Date 10/15/2009	OHED) (USE		



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 08340 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography and found to contain 0.1211 percent (w/vol) ethyl alcohol. The expiration date for this lot number is October 15, 2009 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of **0.10** percent BAC.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204089 06/02/09

resting officer: OWENSZRZD OFFICER I.D.: 813 PERMIT NUMBER: 920086 EXPIRATION DATE: 04/22/11 MISCELLANEOUS DATA:

--- SUPERVISOR MODE ---

SUGARIZ TEST	.000	61:10
BLANK TEST	VERIFIED	01:10
INTERNAL STANDARD	.098	01:10
EXTERNAL STANDARD	.000	01:11
BLANK TEST	.098	01:11
EXTERNAL STANDARD	.000	01:12
BLANK TEST	.098	01:12
EXTERNAL STANDARD	•	01:13
BLANK TEST	.000	6)1.20

N = 3SIM. = .1AVG. = .098 Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204089 06/02/09 00:56

--- DIAGNOSTIC CHECK ---

OKRY COMPUTER: OKAY PROGRAM: HEATERS SAMPLE CHAMBER: 50c OKAY FLOW DETECTOR: PUMP OKAY HIGH SPEED: OKAY DETECTOR: OKAY FILTERS: OKRY QUARTZ STANDARD: OKRY CALIBRATION:

PRINTER TEST !"#\$%&'()*+,-./0123456789:;<=>?@ABCDEFG HIJKLMNOPQRSTUVWXYZE\]^_^abcdefghijklmno pgrstuvwxyz(1)+

or Signature_____ n recycled paper with agri-based inks

CMSU 2208-02

Operator Signature 77/K Printed on recycled paper with agri-based links

CMSU 2208-02

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204089 06/02/09

ARREST TIME: 00:00
SUBJECT NAME:
TEST
DOB: 01/01/01 SEX: M
STATE/D.L.: NA/1234
ARRESTING OFFICER:
N/A
OFFICER I.D.:
TESTING OFFICER:
OWENS/R/D
OFFICER I.D.: 813

OFFICER I.D.: 813
PERMIT NUMBER: 920086
EXPIRATION DATE: 04/22/11
MISCELLANEOUS DATA:

--- BREATH ANALYSIS ---

BLANK TEST

.000

01:17

INTERNAL STANDARD

VERIFIED

01:17

RADIO INTERFERENCE

Operator Signature_

Printed on recycled paper with agri-based inks

CMSU 2208-02

State of Missouri DEPARTMENT OF HEALTH



MO 580-0771 (7-88)

PERMIT TYPE II



Lab. 4 (A7-88)

RICHARD D. OWENS

is hereby authorized to instruct a calibrate, perform field repairs, and	and supervise operators, train instructors, inspect, operate the following breath analyzer(s):		
DATAMASTER for the determination of the alcoholic content of blood from a sample of expired (alveolar) air. Issued under the provisions of sections 577.020 through 577.041, RSMo 1986.			
920086	Director of State Public Health Laboratory		
Expires 04/22/2011	Margart T. Brundly		
	Director, Department of Health		