

MISSOURI DEPARTMENT OF HEALTH STATE PUBLIC HEALTH LABORATORY

LIATH ALCOHOL PROBABL

DATAMASTER MAINTENANCE REPORT

A SHAPE	DATAMASTER MAINT			- R 9" .		
Complete is repaired	this report in duplicate at the . Send copy to Department of	e time of the regu Health; retain origi	lar monthly preventive mainteinal in department file.	nance check, a	nd whenever instrument	
DATAMASTER			(midd diene		inspection 31-2009	
I	INSTRUMENT (STREET AND CITY)			TIME OF	INSPECTION	
CHECKLIS	ST: Place a check (✓) to the	left of each item if t	OATN SHELBYVILLE, found to be satisfactory or if or ust be corrected before using in	erating within	901 HRS established limits. (Write	
	NOSTIC CHECK (PRINTOUT		ast be concered before asing in	ottument.		
	OMPUTER :		DETECTOR			
U PI	ROGRAM		D FILTERS			
HEATERS SAMPLE CHAMBER 49 °C			QUARTZ STANDARD			
FLOW DETECTOR			CALIBRATION			
D PL	JMP HIGH SPEED		PRINTER	PRINTER		
INDIC	ATOR LIGHTS					
TIME /	AND DATE					
☑ SIMUL	ATOR TEMPERATURE (34 °	C ± 0.2°C)				
L CALIB	RATION CHECK -					
spread			tests must be within \pm 5% of to the standard solution being			
0.0	00% STANDARD - MUST REA 40% STANDARD - MUST REA 7 ONE STANDARD IS TO BE	D BETWEEN 0.038	% and 0.042% INCLUSIVE			
TEST 1	, 102	TEST 2	,102	EST3 🏴	, 102	
PERFO	PRM R.F.I. TEST (PRINTOUT	ATTACHED)				
NUMBI FOLLO			REPORT, AND NUMBER OF E	BREATH TESTS	IN EACH RANGE AS	
REFUSALS	6 (004)	(.0509)	(.1014)	.19)	(Over .19)	
and within	w parts and describe any altestablished limits (use other s RUMENT MEETS DIS	ide if necessary)	tion that was made to restore	the instrument	to operate satisfactorily	
			BOTRE# 916 EXP.	10-15-2009		
	F127232			·	,	
INSPECTING SIGNATURE	NG OFFICER		PRINT NAME			
► THE	JJ 173		E. F. BR	JWN .		
TYPE II PERMIT I	NUMBER/EXPIRATION DATE		TELEPHONE NUMBER	85.213	٦	



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 08340 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography and found to contain 0.1211 percent (w/vol) ethyl alcohol. The expiration date for this lot number is October 15, 2009 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcoholanalysis instrument reading of 0.10 percent BAC.

The alcoho and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOUR) STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204098 05/31/09

TESTING OFFICER:
BROWNEVF
OFFICER 1.D.: 173
PERMIT MUMBER: 820056
EXPIRATION DATE: 02/25/10
MISCELLANEOUS DATA:

-- SUPERVISOR MOUE ---

SLAMK TEST	. 066	13:50
INTERNAL STANDARD	VERIF (ED	19:00
EXTORMI STATIANT	. 182	19:50
BLEUF TEST	. 896	19:69
EXTERNAL STAINARD	- 192	19199
BLAW TEST	. 666	
EXTERNAL STATURED	.168	19:18
BLAMK TEST	. POH	4 (15 9 1

M = 3 SIM. = .1 AVG. = .192

Operator Signature TR 22 13

Face This Side Down - This Edge In First

BAC DataMasterEvidence Ticket

MISSOURI STATE HIGHWAY PAIROL BAC DATAMASTER SERIAL HUMBER 204098 05/31/89 19:01

--- DIRBHOSTIC CHECK ---

COMPUTER:	OKAY
PAGGROM:	GGAY
HEATERS SAMPLE CHAMBER:	455
FLOW DETECTOR:	OKAY
PUMP MIGH SPEED:	CKAY
DETECTOR:	OKRY
FILTERS	CERY
OUARTZ STANDARN	OKHY
) CHLIBATIONS	ONAY

PRINTER TEST

!"##X&'()##.-./0123456709*;(=>?@HBCDEFG HTUKLNMOPORSTUVW(YZENIO_*abcdefghi.kimno prinstuvwxyz(I)*"

Operator Signature TR 6 2 5 5 5

2208-0

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PRIROL BAC DATAMASTER SERIAL MUMBER 204898 95731709

ARREST TIME: 19:00 SUBJECT HAME: TEST DOB: 01/01/01 SEX: M STATE/D.L.: MO/TEST ARRESTIMG OFFICER; TEST OFFICER I.D.: 173 TESTIMO OFFICERS BROWNZEZF OFFICER L.D.: 173 PERMIT HUMBER: 820056 EXPIRATION DATE: 02/25/10 MISCELLAMEOUS DATE:

--- BREATH AWALYSIS ---

RADIO INTERPEREDCE

Operator Signature TH. Carte 138

2208-02

State of Missouri DEPARTMENT OF HEALTH



TYPE II



E.F. BROWN

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field repairs, and operate the following breath analyzer(s):

DATAMASTER

for the determination of the alcoholic content of blood from a sample of expired (alveolar) air. Issued under the provisions of sections 577.020 through 577.041, RSMo 1986.

02/25/08	Time C. Ostant			
Number 820056	Director of State Public Health Laboratory			
02/25/2010	District Additional Control of Manufacture and C			
•	Director, Department of Health			

MO 680-0771 (7-88)