

MISSOURI DEPARTMENT OF HEALTH STATE PUBLIC HEALTH LABORATORY

CORATH ALCOHOL PRODUCAL

DATAMASTER MAINTENANCE REPORT

Complete this report in duplicate at the time of the regular is repaired. Send copy to Department of Health; retain original	monthly preventive maintenance coin department file.	heck, and whenever instrumen
DATAMASTER SN 204105	S 400	DATE OF INSPECTION 05-31-09
LOCATION OF INSTRUMENT (STREET AND CITY)		TIME OF INSPECTION
SULLIVAN COUNTY SHERIFF'S O	PFICE	1031
CHECKLIST: Place a check (<) to the left of each item if four in observed values where determined.) Unchecked items must lead to the control of the control	nd to be satisfactory or if operating be corrected before using instrumer	within established limits. (Write at.
DIAGNOSTIC CHECK (PRINTOUT ATTACHED)		
COMPUTER .	DETECTOR	
PROGRAM	FILTERS	
HEATERS SAMPLE CHAMBER 49 °C	☐ QUARTZ STANDARD	
FLOW DETECTOR	☐ CALIBRATION	
PUMP HIGH SPEED	PRINTER	
M INDICATOR LIGHTS		
TIME AND DATE 1031 05-31-09		
SIMULATOR TEMPERATURE (34 °C ± 0.2°C) 34,0°	<i>(</i> .	
CALIBRATION CHECK -		
Run three tests using a standard solution. All three test	s must be within \pm 5% of the sta	ndard value and must have a
spread of .005 or less. Check the box corresponding to the	ne standard solution being used. (F	PRINTOUT ATTACHED) (USE
RECIRCULATION PUMP)	. 10 405% [NO. 110]	
☑ 0.100% STANDARD - MUST READ BETWEEN 0.095% at ☐ 0.040% STANDARD - MUST READ BETWEEN 0.038% at	nd 0.105% INCLUSIVE nd 0.042% INCLUSIVE	•
(ONLY ONE STANDARD IS TO BE USED PER MAINTENA		
TEST 1 - 102 TEST 2 -	102 TEST 3	,102
PERFORM R.F.I. TEST (PRINTOUT ATTACHED)		•
NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH RANGE AS FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS)		
REFUSALS / (004) / (.0509) /	(.1014) 3 (.1519)	(Over .19)
List any new parts and describe any alteration or modification	that was made to restore the instr	ument to operate satisfactorily
and within established limits (use other side if necessary) GOTH REPCO SIMULATOR MODEL	SUMAC SECURI & III	19
401 # 08340	- 370/C- XXINC - 11-	<u> </u>
EXP. 10-15-09		· · · · · · · · · · · · · · · · · · ·
100% VAPOR		
INSTRUMENT IS OPERATING WI	THIN DHSS SPECIA	FICATIONS
INSPECTING OFFICER	PRINT NAME	
El B. S. Ream	CPL.B.D. REAM	
TYPE II PERMIT NUMBER/EXPIRATION DATE	TELEPHONE NUMBER	
820065 02-25-10	(660) 385-2132	



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 08340 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography and found to contain 0.1211 percent (w/vol) ethyl alcohol. The expiration date for this lot number is October 15, 2009 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- 2°C, this solution will give a breath alcoholanalysis instrument reading of 0.10 percent BAC.

The alcoho and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

Face This Side Down - This Edge In First

BAC DataMasterEvidence Ticket

MISSOURI STATE HIGHWAY PATROL EAC DATAMASTER SERIAL NUMBER 204105 05/31/09 18:31

--- DIAGNOSTIC CHECK ---

CUMPUTERS

OKAY

FRUGRAM

(REY

HERIERS

SAMPLE CHAMBERS

g Şu-

FLUM DETECTORS

OKAY

PHAR

HIGH SPEED:

OKHY

DETECTORS

OKHY

FILTERS

(1)(1)(1)

OURREZ STRABARIO

THEY

CALIBRATION:

OKAY

PRINTER TEST

!"##%%:(()**,-./0123456789:;(=)?@ABCUCF6 HI.KLMNOP@RSTUVWWYZE\I^_^abcdefghijk}mno pqrstuvwwgz{I}+~

Operator Signature_

on recycled paper with agri-based inks

ONOU 0000 00

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204105 05/31/09

--- BREATH AMALYSIS ---

RADIO INTERFERENCE

Face This Side Down - This Edge In First

BAC DataMasterEvidence Ticket

MISSOURI STATE HIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204105 05/31/09

TESTING OFFICER:
REAM/B/D
OFFICER I.D.: 356
PERMIT NUMBER: 829065
EXPIRATION DATE: 02/85/10
MISCELLAMEOUS DATA:

--- SUPERVISOR NODE ---

RLANK TEST	, 900	19637
INTERNAL STRUCKS	VER (F) ED	13:07
EXTERNAL STANDARD	. 162	18439
BLANK TEST	, Geni	.10:39
EXTERME STANDARD	. T. (1)	19141
BLAPK TEST	. due	19:41
EXTERNAL STANDARD	.102	10:43
BLARK TEST	. 565	10143

n = 3 Sim. = .1 Tevs. = .102

Operator Signature 8 8 8 Rocers # 350 CM

Printed on recycled paper with agri-based inks

Operator Signature

Printed on recycled paper with agri-based inks

CMSH 2208-02

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

MISSOURI STATE HIGHWAY PATROL.: BAC DATAMASTER SERIAL HUMBER 204105 05/31/09

ARREST TIME: 10:15
SUBJECT MAME:
DOE/JOHN/E
DOB: 01/01/70 SEX: M
STATE/D.L.: MO/999999999
ARRESTIME OFFICEF:
REAM/R/D
OFFICER 1.D.: 356
TESTIME OFFICER:
REAM/B/D
OFFICEP 1.D.: 356
PERMIT NUMBER: 820065
EXPIRATION DATE: 02/25/10
MISCELLAMEOUS DATA:

.... DREATH AMALYSIS ---

REDIO INTERFERENCE

Face This Side Down - This Edge In First

BAC DataMasterEvidence Ticket

MISSOURI STATE MIGHWAY PATROL BAC DATAMASTER SERIAL NUMBER 204105 05/31/09

TESTING OFFICER:
REMMUBUD
OFFICER I.D.: 356
PERMIT NUMBER: 820065
EMPIRATION DATE: 02/25/10
MISCELLAMEOUS DATA:

--- SUPERVISOR NODE ---

ELMIK TEST	. (984)	19132
INTERNAL STRABBED	VER (F) (I)	18:37
EXTERMAL STAMOARD	.168	19:59
BLANK TEST	, 1961 <u>i</u> 9	19:30
EXTERNAL STAMARD	. 195	19141
PLEAR TEST	. 609	19:41
EXTERNAL STANDARD	100	10:45
THANK TEST	. 19(36)	18849

n = 3 Sin. = .1 My6, = .102

Operator Signature 68 B. D. Roccom #3

Printed on recycled paper with agri-based inks

Operator Signature

Printed on recycled paper with agri-based inks

CMSH 2209-0

State of Missouri DEPARTMENT OF HEALTH



TYPE II



BRADLEY REAM

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field repairs, and operate the following breath analyzer(s):

DATAMASTER

for the determination of the alcoholic content of blood from a sample of expired (alveolar) air. Issued under the provisions of sections 577.020 through 577.041, RSMo 1986.

02/25/08	Tui C. Oolanh
820065	Director of State Public Health Laboratory
Expires 02/25/2010	To the state of th
	Director, Department of Health

MO 680-0771 (7-88)