



MISSOURI DEPARTMENT OF HEALTH
STATE PUBLIC HEALTH LABORATORY
DATAMASTER MAINTENANCE REPORT

BREATH ALCOHOL PROGRAM

JUN-2 2009

Complete this report in duplicate at the time of the regular monthly preventive maintenance check, and whenever instrument is repaired. Send copy to Department of Health; retain original in department file.

DATAMASTER SN 940444	DATE OF INSPECTION 04-27-09
LOCATION OF INSTRUMENT (STREET AND CITY) 1200 Main St Grandview Mo 64030	TIME OF INSPECTION 12 09 Hours

CHECKLIST: Place a check (✓) to the left of each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unchecked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	
<input checked="" type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER 49 °C	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER
<input checked="" type="checkbox"/> INDICATOR LIGHTS	
<input checked="" type="checkbox"/> TIME AND DATE	
<input checked="" type="checkbox"/> SIMULATOR TEMPERATURE (34 °C ± 0.2°C) 34.0C	
<input type="checkbox"/> CALIBRATION CHECK - Run three tests using a standard solution. All three tests must be within ± 5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) (USE RECIRCULATION PUMP) <input type="checkbox"/> 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE <input type="checkbox"/> 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)	

TEST 1 <input checked="" type="checkbox"/> .098	TEST 2 <input checked="" type="checkbox"/> .098	TEST 3 <input checked="" type="checkbox"/> .098
--	--	--

<input checked="" type="checkbox"/> PERFORM R.F.I. TEST (PRINTOUT ATTACHED)						
<input checked="" type="checkbox"/> NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH RANGE AS FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS)						
REFUSALS 12	(0-.04) 3	(.05-.09) 3	(.10-.14) 1	(.15-.19) 7	(Over .19) 1	

List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary)

Instrument Meets D.O.H. Specs.
Used Certified .10% solution from repco Marketing
Lot # 08002 Manufacture Date 10-14-08 Exp Date 10-13-10

INSPECTING OFFICER	
SIGNATURE 	PRINT NAME Ricky L. VanSickie
TYPE PERMIT NUMBER/EXPIRATION DATE 820191 06-09-10	TELEPHONE NUMBER 816-316-4900

CERTIFICATE OF ANALYSIS

Random samples of lot number 08002 of Alcohol Certified Solution for simulator were analyzed by an independent laboratory utilizing a gas chromatograph and found to contain .1209 gms/dl wt. /vol. ethyl alcohol.

The alcohol and distilled water used in the solution were found to be free of any interfering substance.

When used in a calibrated simulator, operating at 34 c +/- .2 c, this solution will give an alcohol breath test instrument reading of .100 percent BAC +/- 2% or .002 BAC (whichever is greater).

The expiration date for this lot number is October 13, 2010 at 11:59PM.

This document is a true representation of the original Certificate of Analysis.



Cecil B. Garner, President
RepCo Marketing, Inc.

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI
GRANDVIEW POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 940444
05/27/09
12:16

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY

PROGRAM (04-07-2009): OKAY

HEATERS

SAMPLE CHAMBER: 49c

FLOW DETECTOR: OKAY

PUMP

HIGH SPEED: OKAY

DETECTOR: OKAY

FILTERS: OKAY

QUARTZ STANDARD: OKAY

CALIBRATION: OKAY

PRINTER TEST

!"#\$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJ
KLMNOPQRSTUVWXYZ[\]^_`abcdefgijklmnop
qrstuvwxyz{|}~

OPERATOR SIGNATURE

CARD STK #
60036

REORDER ALL SUPPLIES FROM N.P.A.S.
2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (NPAS)

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI
GRANDVIEW POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 940444
05/27/09

TESTING OFFICER:

VANSICKLE/RICKY/L
OFFICER I.D.: 0168
PERMIT NUMBER: 820191
EXPIRATION DATE: 06/09/10
MISCELLANEOUS DATA:

--- SUPERVISOR MODE ---

BLANK TEST	.000	12:09
INTERNAL STANDARD	VERIFIED	12:09
EXTERNAL STANDARD	.098	12:10
BLANK TEST	.000	12:10
EXTERNAL STANDARD	.098	12:11
BLANK TEST	.000	12:11
EXTERNAL STANDARD	.098	12:12
BLANK TEST	.000	12:13

N = 3
SIM. = .1
AVG. = .098

OPERATOR SIGNATURE

CARD STK #
60036

REORDER ALL SUPPLIES FROM N.P.A.S.
2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (NPAS)

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI
GRANDVIEW POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 940444
05/27/09

ARREST TIME: 12:00

SUBJECT NAME:

TEST

DOB: 11/22/33 SEX: M

STATE/D.L.: MO/123456789

ARRESTING OFFICER:

VANSICKLE/RICKY/L

OFFICER I.D.: 0168

TESTING OFFICER:

VANSICKLE/RICKY/L

OFFICER I.D.: 0168

PERMIT NUMBER: 820191

EXPIRATION DATE: 06/09/10

MISCELLANEOUS DATA:

--- BREATH ANALYSIS ---

BLANK TEST	.000	12:24
INTERNAL STANDARD	VERIFIED	12:25
RADIO INTERFERENCE		

OPERATOR SIGNATURE

CARD STK #
60036

REORDER ALL SUPPLIES FROM N.P.A.S.
2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (NPAS)