

STATE PUBLIC HEALTH LABORATORY

DATAMASTER MAINTENANCE REPORT

complete this report in duplicate at the time of repaired. Send copy to Department of Health;	of the regular mo	nthly preventive mainte		1	
repaired. Send copy to Department of Health,	retain onginarii C		DATE OF INS	COHOL PROS	
950123 La Monte Police				23/2009 PECTION 1 -8 7	
304 West Front Street.La.Mont	e	- taleston or if o	1/56	abiliahad limits (Write)	
CHECKLIST: Place a check (<) to the left of elements of the left of elements of the control of t	ach item it found ted items must be	corrected before using I	nstrument	Total Control of the	
DIAGNOSTIC CHECK (PRINTOUT ATTAC		Abelle,	<u>,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,</u>	114.11	
COMPUTER		DETECTOR			
PROGRAM		FILTERS			
HEATERS SAMPLE CHAMBER 49	<u>L.c</u>	QUARTZ STANDAR	ID.		
FLOW DETECTOR		CALIBRATION	المتعادلة في إنها إن إنجاب مساعد		
PUMP HIGH SPEED		PRINTER			
INDICATOR LIGHTS					
TIME AND DATE					
SIMULATOR TEMPERATURE (34 °C ± 0.2	200) 34°C				
Run three tests using a standard solution spread of .005 or tess. Check the box connected the connected purple of .100% STANDARD - MUST READ BETONLY ONE STANDARD IS TO BE USED	rresponding to th TWEEN 0.095% an TWEEN 0.038% an	e atendard solution beind 0.105% INCLUSIVE Id 0.042% INCLUSIVE	of the standard va ig used. (PRINTOL	lue and must have a IT ATTACHED) (USE	
	EST 2	096	TEST 3 .096		
PERFORM R.F.I. TEST (PRINTOUT ATTA	CHED)				
NUMBER OF REFUSALS, SINCE LAST MATERIAL FOLLOWS: (DO NOT INCLUDE SIMULA		PORT, AND NUMBER OF	F BREATH TESTS	N EACH RANGE AS	
REFUSALS 004) 0 (.05-	(00,-	(.1014)	1619)	(Over .19).	
List any new parte and describe any atteration and within established limits (use other side if the content of		that was made to resto	re the instrument to	operate satisfactorily	
				• . •	
SUPPLIER: Guth Labs	CONCENTRA'	rion: 0.10 EXPIRE	81 _{8/11/2009} LOT	08280	
INSPECTING OFFICER		PRINT NAME			
► (Journ Cll	·	TELEPHONE NUMBER	James E. Uln	<u>n</u>	
TYPE II PERMIT NUMBER/EXPIRATION DATE 720199 / 09/27	//2009	1	660-347-5106		
	AN ADULT OPPOSTIBATIVA	CHIMATON ACTION PUPLOYER		(Jab11)	



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 08280 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography and found to contain 0.1212 percent (w/vol) ethyl alcohol. The expiration date for this lot number is August 11, 2009 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.10 percent BAC.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

CAMONTE POLICE DEPARTMENT

SAC IATAMASTER SERIAL NUMBER 950123 06/03/09

OFFICER I.D.: 781 PRESET NUMBER: 700199 THISTING OFFICER: NSCELLANEOUS DATA: EXPERATION DATE: 89/25/89

SCHEREISON MODE ---

20年2月 - 1987年 - 1987	ļā g	语声品
 2008 2008 2008	 	. 988 . 986 . 996 . 996
20 20 47 20 20 41 10 10 10 14 14 14	40	

60 # # 60 # # 60

Face This Side Down - This Edge In First

- ;-

BAC DataMaster Evidence Ticket

LAMBUTE POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 950123

OFFICER I.D.: 781 PERMIT NUMBER: 720199 EXPLARTION DATE: 09/25/09 RISCELLANEOUS DATA: TESTING OFFICERS DOB: 11/11/11 SEX: P STATE/D.L.: MD/123456789 ARRESTING OFFICER: ERREST TIME: 18:88 产艺生的

SECTE PERSON

REDIC LATERIERONDE

. 2000 : 8

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

CHMONTE POLICE DEPARTMENT:

SAC DATAMASTER SERIAL NUMBER 950123

DIAGNOSTIC CHECK ---

STRREAM STREETS	THE PROOF	DETECTOR:	HIGH SPEED: PURD	FLOW DETECTOR:	CHAPLE CLARGER:	PROBRAM	COMPUTER:
OKHY	VEXIO	DKBY .	OKRY	OKAY	490 0	DK9Y	9 , 97

Persetucionem(1)中本 #117C1#50PQRSTUV8XYZ[737]", abcdefeni3k1mno CHICKATION

9499

Printed on recycled paper with agre-based inks Operator Signature

CMSU 2208-02

Printed on recycled paper with Operator Signature

CMSU 2208-02

Printed on recycled paper with agri-based into

CMSU 2208-02

Operator Signature