

MISSOURI DEPARTMENT OF HEALTH STATE PUBLIC HEALTH LABORATORY DATAMASTER MAINTENANCE REPORT

BREATH ALCOHOL PROGRAM

JUN - 9 7009

48(4)	• •	001	v ~ 3	ZIIII4	
Complete this report in duplicate at the time of the regula is repaired. Send copy to Department of Health; retain original	r monthly preventive ma al in department file.	Intenance ch ∪i135 \$1,	eck an	nd whonover in	strument
DATAMASTER SN 201300			DATE OF	INSPECTION	
LOCATION OF INSTRUMENT (STREET AND CITY)				NSPECTION	
200 N. Bourbeuse St. James, Mo. 65559 CHECKLIST: Place a check (<) to the left of each item if fo	und to be satisfactory or	if operating v	1416hr vithio e		s (Write
in observed values where determined.) Unchecked items mus	t be corrected before usin	g instrument			
☑ DIAGNOSTIC CHECK (PRINTOUT ATTACHED)					
☑ COMPUTER	☑ DETECTOR				
☑ PROGRAM	☑ FILTERS				
☑ HEATERS SAMPLE CHAMBER 49 °C	☑ QUARTZ STANDARD				
☑ FLOW DETECTOR	☑ CALIBRATION				
☑ PUMP HIGH SPEED	PRINTER				
☑ INDICATOR LIGHTS					
☑ TIME AND DATE		en e			
SIMULATOR TEMPERATURE (34 °C ± 0.2°C) 34.0 deg	rees C				****
☑ CALIBRATION CHECK -					
Run three tests using a standard solution. All three te	sts must be within \pm 5%	of the stan	dard va	alue and must	have a
spread of .005 or less. Check the box corresponding to RECIRCULATION PUMP)	the standard solution be	ing used. (Pl	JOTAIF	JT ATTACHED) (USE
0.100% STANDARD - MUST READ BETWEEN 0.095%	and 0.105% INCLUSIVE				
☐ 0.040% STANDARD - MUST READ BETWEEN 0.038% (ONLY ONE STANDARD IS TO BE USED PER MAINTEN					
TEST 1 .101 TEST 2	.101	TEST 3		.101	
PERFORM R.F.I. TEST (PRINTOUT ATTACHED)					
☑ NUMBER OF REFUSALS, SINCE LAST MAINTENANCE F	REPORT AND MILIABED (NE DDEATH 1	reere	IN EACH DANK	
FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS)	ici Offi, AND NONBEST C	A DULTIL	COTO	IN EACH HANC	IE NO
REFUSALS 1 (004) 0 (.0509) 0	(.1014) 0	(.1519)	0	(Over .19)	0
List any new parts and describe any alteration or modification and within established limits (use other side if necessary)	on that was made to resto	ore the instru	ment to	o operate satisf	actorily
Instrument operating within DOH specifications.					
		V-V-A-7-100-100-10-10-10-10-10-10-10-10-10-10-1			
REPKO MKTG .100% solution Exp. October 13, 2010 Lo	ot#08002			-	
NSPECTING OFFICER	Service and California et ale and de				
IGNATURE C: M. C.	PRINT NAME Timothy C. McDona	14			
YPE II PERMIT NUMBER/EXPIRATION DATE (20160 05/12/2010	TELEPHONE NUMBER	iu			
20100 03/12/2010	(573) 265-7012				

CERTIFICATE OF ANALYSIS

Random samples of lot number <u>08002</u> of Alcohol Certified Solution for simulator were analyzed by an independent laboratory utilizing a gas chromatograph and found to contain <u>.1209</u> gms/dl wt. /vol. ethyl alcohol.

The alcohol and distilled water used in the solution were found to be free of any interfering substance.

When used in a calibrated simulator, operating at 34 c \pm -.2 c, this solution will give an alcohol breath test instrument reading of .100 percent BAC \pm -.2% or .002 BAC (whichever is greater).

The expiration date for this lot number is October 13, 2010 at 11:59PM.

This document is a true representation of the original Certificate of Analysis.

Cecil B. Garner, President

RepCo Marketing, Inc.

FACE THIS SIDE DOWN - THIS EDGE IN FIRS

<u>าคา เพลง ซาการสุดสุดสุดสุดสุดสุดสุดสุดสุด</u>

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI ST. JAMES POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201300

ARREST TIME: 13:00 SUBJECT NAME:

MCLOVIH

DOB: 89/01/90

SEM: M

STATE/D.L.: MI/1234567809

ARRESTING OFFICER:

MCDONALDZTIMOTHYZČ

DEFICER I.D. 9 1192

TESTING OFFICER:

MCDOWALD/TIMOTHY/C

OFFICER I.D. = 1102

PERMIT NUMBER: 820160

EXPIRATION DATE: 05/12/10

MISCELLAMEOUS DATA:

REI TEST

--- BREATH AMALYSIS -

EADIO INTERFERENCE

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI ST. JAMES POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201300

TESTING OFFICER:

MCDOMALI/TIMOTHY/C
OFFICER I.D.: 1102
PERMIT MUMBER: 820160
EXPIRATION DATE: 05/12/10
MISCELLAMEOUS DATA:

--- SUPERVISOR MODE ---

BLANK TEST	. 665	14:16
INTERMAL STAMBARD	VERIFIED	14:16
EXTERMAL STAMBARD	. 181	14116.
BLANK TEST	. Gliği	14:17
EXTERNAL STAMDARD	.101	1 4 4 1 1
BLANK TEST	. G68	14:18
EXTERNAL STANDARD	. 101	14118
BLANK TEST	. 688	14:18

N = 3 SIN. = .1 AVG. = .101 _

OPERATOR SIGNATURE

CARD STK # 60036

REORDER ALL SUPPLIES FROM N.P.A.S.

2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (N

OPERATOR SIGNATURE,

CARD STK # 60036

REORDER ALL SUPPLIES FROM N.P.A.S.

2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (NPAS)

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI ST. JAMES POLICE DEPARTMENT

PAC DATAMASTER SERIAL NUMBER 201300 06/05/09 14:20

--- DIAGNOSTIC CHECK ---

COMPUTERS

OKAY

PROGRAM:

ÜKAY

HEATERS

SAMPLE CHAMBER:

495

FLOW DETECTOR:

OKAY

41191

HIGH SPEED:

OKAY

DETECTOR:

OKAY

FILTERS:

OKAY

QUARTZ STANDARD:

OKAY

CALIEMATION:

OKAY

PRINCER TEST

!"#\$%&(()*+,-./0123456789:;<=>?@ABCDEFG HIUKLMNOPQRSTUVWXYZ[\1^__'abcdefghijklmno pqrstuvwxyz(1)="

OPERATOR SIGNATURE /

CARD STK # 60036

REORDER ALL SUPPLIES FROM N.P.A.S.

2260 NORTH MAIN, MANSFIELD, OH 44903 419-526-6727 (NPAS)

State of Missouri DEPARTMENT OF HEALTH



PERMIT TYPE II



TIMOTH	HY MCDONALD
is hereby authorized to instruct calibrate, perform field repairs, and	and supervise operators, train instructors, inspect, i operate the following breath analyzer(s):
DATAMA	ASTER
for the determination of the alcohol air. Issued under the provisions of s	ic content of blood from a sample of expired (alveolar) sections 577.020 through 577.041, RSMo 1986.
Dale05/12/08	Time C. Dolland
Number 820160	Director of State Public Health Laboratory
Expires 05/12/2010	
MO 680-0771 (7-88)	Director, Department of Health Lab. 4 (R7-88)
	L40. 4 (nr-80)