



MISSOURI DEPARTMENT OF HEALTH
STATE PUBLIC HEALTH LABORATORY

DATAMASTER MAINTENANCE REPORT

BREATH ALCOHOL PROGRAM

JUN - 8 2009

Complete this report in duplicate at the time of the regular monthly preventive maintenance check, and whenever instrument is repaired. Send copy to Department of Health; retain original in department file.

DATAMASTER SN		DATE OF INSPECTION
BAC Datamaster 95 0089		06/07/09
LOCATION OF INSTRUMENT (STREET AND CITY)		TIME OF INSPECTION
104 W. 3rd street Lawson Mo. 64062		1422

CHECKLIST: Place a check (✓) to the left of each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unchecked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	
<input checked="" type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER 50 °C	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER
<input checked="" type="checkbox"/> INDICATOR LIGHTS	
<input checked="" type="checkbox"/> TIME AND DATE	
<input checked="" type="checkbox"/> SIMULATOR TEMPERATURE (34 °C ± 0.2°C) (34.1 °C)	

CALIBRATION CHECK -

Run three tests using a standard solution. All three tests must be within \pm 5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) (USE RECIRCULATION PUMP)

- 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE
 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE
 (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

TEST 1	0.102%	TEST 2	0.102%	TEST 3	0.103%
--------	--------	--------	--------	--------	--------

PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

NUMBER OF REFUSALS, SINCE LAST MAINTENANCE REPORT, AND NUMBER OF BREATH TESTS IN EACH RANGE AS FOLLOWS: (DO NOT INCLUDE SIMULATOR TESTS)

REFUSALS	1	(0-.04)	0	(.05-.09)	0	(.10-.14)	0	(.15-.19)	0	(Over .19)	0
----------	---	---------	---	-----------	---	-----------	---	-----------	---	------------	---

List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary)

Instrument conforms to DOH-SS STANDARDS

MANUFACTURER: PEPCO INSPECTING OFFICER:	CONCENTRATION: 0.100% EXPIRES: 06/11/09 LOT #: 08001
SIGNATURE 	PRINT NAME P. Brian Lefever
TYPE II PERMIT NUMBER/EXPIRATION DATE 820 221 07/23/10	TELEPHONE NUMBER Lawson P.D. (816) 580-7210

RepCo MARKETING INC.

3101-180 STONYBROOK DRIVE
RALEIGH, NC 27604
919-876-5480

CERTIFICATE OF ANALYSIS

MANUFACTURER AND SUPPLIER: RepCo Marketing, Inc.

LOT NUMBER: 08001

EXPIRATION DATE: August 11, 2009 at 11:59 p.m.

RepCo Marketing, Inc. certifies the following:

RepCo Marketing, Inc. manufactured, tested and supplied Lot Number 08001 of Alcohol Certified Solution for simulators. Random samples of said lot number were analyzed utilizing a gas chromatograph and found to contain 0.1216 gms/dl wt/vol ethyl alcohol in aqueous solution (i.e. ethanol).

The alcohol and distilled water used in the solution were found to be free of any interfering substance.

This solution will produce a vapor alcohol value of 0.100% +/-3% when heated to 34 Degrees Celsius +/-0.2 Degrees Celsius in a simulator.

The date of manufacture for this lot number is February 12, 2008.
The expiration date for this lot number is August 11, 2009 at 11:59 p.m.

This document is a true representation of the original Certificate of Analysis.

Cecil B. Garner

Cecil B. Garner, President
RepCo Marketing, Inc.

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI
LAWSON POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 9500089
06/03/09
14:22

— DIAGNOSTIC CHECK —

COMPUTER: OKAY

PROGRAM: OKAY

HEATERS: SAMPLE CHAMBER: SEC

TESTING OFFICER:

LAFAVOR

OFFICER I.D.: 200

PERMIT NUMBER: 8200221

EXPIRATION DATE: 07/23/10

MISCELLANEOUS DATA:

— SUPERVISOR MODE —

BLANK TEST	000	14:25
INTERNAL STANDARD	VERIFIED	14:25
EXTERNAL STANDARD	.102	14:25
BLANK TEST	.000	14:26
EXTERNAL STANDARD	.102	14:26
BLANK TEST	.000	14:27
EXTERNAL STANDARD	.102	14:27
BLANK TEST	.000	14:28

N = 3
STM. = .1
RVS. = .1023

PRINTER TEST
!#%& '()@12345678@:;<=^@HBCDEFS
HIJKLMNOPQRSTUVWXYZ\j`_abdefhijklmn
PQRSTUVWXYZ\j`_abdefhijklmn

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster

Evidence Ticket

STATE OF MISSOURI
LAWSON POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 9500089
06/03/09

TESTING OFFICER:

LAFAVOR

OFFICER I.D.: 200

PERMIT NUMBER: 8200221

EXPIRATION DATE: 07/23/10

MISCELLANEOUS DATA:

— BREATH ANALYSIS —

BLANK TEST	000	14:31
INTERNAL STANDARD	VERIFIED	14:31
RADIO INTERFERENCE		

OPERATOR SIGNATURE

J. S. L.

CARD STK # REORDER ALL SUPPLIES FROM N.P.A.S.
60036 P.O. BOX 1435, MANSFIELD, OH 44901

OPERATOR SIGNATURE

J. S. L.

CARD STK # REORDER ALL SUPPLIES FROM N.P.A.S.
60036 P.O. BOX 495, MANSFIELD, OH 44901

State of Missouri
DEPARTMENT OF HEALTH



P E R M I T
TYPE II



P. BRIAN LAFAYOR

is hereby authorized to instruct and supervise operators, train initiators, inspect, calibrate, perform field repairs, and operate the following breath analyzer(s):

DATAMASTER

for the determination of the alcoholic content of blood from a sample of expired (alveolar) air, issued under the provisions of sections 677.020 through 677.041, 119 Mo. 1980.

Date 07/23/08
Number 820221
Expires 07/23/2010
Initials PL
Mo. Reg. No. 07717709

Eric C. Oberle
Director of Ohio Public Health Laboratory
 Director, Department of Health
 Lab. A (10-00)